Sample Preparation

for the Laser Particle Sizer ANALYSETTE 22

Sample preparation

1. Sample division

2. Sample preparation

2.1 Dispersion

2.2 Tips and tricks

3. Materials and suitable measuring liquids

Sample preparation

It is always astounding that the market for highly precise, fully automatic analytical instruments is permanently growing, while no emphasis is placed on the equally important sample preparation or sampling.

One of the most common mistakes during analysis is already made during the sampling stage and drawn into the analytical procedure right from the beginning. Again and again, it surprises process technicians how often samples are drawn carelessly and how the achieved analysis results are readily accepted.

Initially the capabilities of the measuring instrument are often doubted when repeat measurements deviate, but the source of the error is deeper though: in an inhomogeneous sampling.

Therefore the results are only reproducible if the analysed sample **represents** the goods to be tested with a high degree of exactness, i.e. the sample taken can be equated with the entire batch.

1. Sample division

The sample for the **Laser Particle Sizer ANALYSETTE 22** (approximately 200 mg – 1 g) shall correspond in the contained particle type and particle distribution with the entire batch of material.

The **Rotary Cone Sample Divider LABORETTE 27** is very well suited for the division of dry laboratory samples or suspensions, since different dividing heads with different division ratios can be selected. Depending on the division head up to 3000 dividing steps per minute can be achieved.

Possible division: 1:8, 1:10 and 1:30

Example: a laboratory sample of 50 ml (g) respectively 100 ml (g), could initially be divided with the dividing head 1:10 and then with the dividing head 1:30, so that for the Laser Particle Sizer ANALYSETTE 22 required amount of 200-400 mg is available.

Now the **representative** sample will be deglomerated.

2. Sample preparation

Via pre-tests it must be determined in which manner the sample material can be moistened and dispersed.

The liquid should possibly completely and spontaneously wet the solid matter.

The additional ultrasonic support (if possible with maximum power) in general reduces the dispersion duration.

- The particles of the solid matter shall be individually and free of agglomerates in the suspension.
- The dispersion condition must be *stable* during the entire measurement, coagulating / flocculating of the particles may not occur.
- Coagulating in a suspension can be recognized by a slight, slushy bottom deposit, which will swim in the sample glass at slight movements of the glass, as a second phase respectively layer and displays a *cloudlike* shape.
- Floating of the sample material on the surface of the liquid is a sign of non-wetted sample material, i.e. the already added share of dispersants like tensides/ wetting agents or salts is too low.
- The sample may not be destroyed during the dispersion, respectively comminuted. With thin, platelet shaped material like mica, kaolin, clay and inorganic salts this is especially important.

Here it is also recommended to paste a small sample amount on a specimen holder with a little bit of liquid and wetting agent. Now under the microscope the particle spectrum and the maximum particle size can be determined.

After the finished dispersion it is also microscopically checked if the coarse particles are still present or were destroyed.

- The vehicle liquid/ measuring liquid must in all cases have a smaller, max. equal specific weight (density) than the solid material to be measured.
- The sample may not partially dissolve, dissolve or swell!

2.1 Dispersion

Unproblematic samples, which submerge without any great effort directly into the surface of the water and possess no large amount of fine share, are added to the dispersion unit as a *solid material* with a spatula portion by portion and after a brief ultrasonic treatment/ dispersion (30-60s) measured reproducibly.

Shows the conducted *double measurement* too large differences in the particle size distribution it may be attributed to several reasons:

- pump speed too low – the coarse material deposits
- the agitator speed is too great – air bubbles are created, respectively air is stirred in
- dispersion duration is too brief

- Increasing fine share
 - Longer ultrasonic treatment necessary
 - Addition of dispersion aids
 - Fine particles stick to the measuring cell glass
 - Clogging due to coarse particle >2 mm
 - Coarse material deposits: specific weight too high
- Seemingly coarser appearing curve
 - Sample swells or flocculates
 - Sample agglomerates
 - Sample is magnetic
- Declining coarse- and fine range
 - Sample dissolves
 - Beam absorption drops

Additionally faulty measurements due to jammed or bent hose connections can occur.

Samples hard to disperse may display the following characteristics: static charge (for example plastics), adhesive forces or cohesive powers – samples tend to conglutinate/ agglutinate – (for example clays, soil samples, kaoline), magnetism, hydrophobic characteristics – water repellent molecule components (for example drugs, medications, toner, graphite, titanium dioxide, waxes), coagulation (for example clays, kaoline, chalk, gypsum).

a) Static charge respectively also hydrophobic characteristics:

Here a spatula amount of the material should be added to a small 50 ml Erlenmeyer flask and then at first 1 (up to 2) drops of a wetting agent (surfactant or diluted surfactant solution) added, then mixed into a paste until the sample is completely wetted. Water is added drop by drop and stirred. The now created suspension (approximately 20 - 30 ml) is dispersed in an ultrasonic bath.

If the sample is already inside the dispersion unit and floats on top of the surface, the sample can be wetted as follows:

With a glass rod or with the tip of a spatula, a *small* drop of the wetting agent (for example Dusazin 901, Teepol, Tween 80 or dish soap) is added/ touched on the surface of the liquid and distributed. Immediately it can be seen, that the formed skin on the surface breaks open and the fine portion enters the suspension.

b) Adhesive powers:

are *clinging* powers of the particle. A reduction of these surface powers can be obtained by creating on the boundary phase areas from solid/ liquid, adsorption layers from surface-active agents or macromolecules. It is considered a covering, a guarding or masking of the solid material respectively wetting. For this reason in most cases for example tetra-Na-diphosphat (sodium pyrophosphate: Na₄P₂O₇) or poly sodium approximately at 0.5-1% are utilized.

c) *Cohesive powers and magnetism:*

The effects of the forces of attraction between atoms or molecules of a body are described as cohesive power or polarity respectively magnetic characteristics. The magnetic characteristics are difficult to eliminate: with only slight magnetism a highly viscous liquid like for example ethylene glycol or a glycerine/ water mixture may be used or the sample is heated above several 100°C. This is hardly possible in a laboratory and seldom realisable. Therefore are such materials with a high degree of magnetism are not suitable for the particle size analysis.

d) *Coagulation:*

is the flocculation of a sample caused by the agglomeration of colloid particles in a suspension. This can happen with too great amounts of solids or an unfavourable pH-range. By adding several drops of undiluted acid (for example hydrochloric acid) prior to adding the sample!! for the acidic range or through diluted lye (for example caustic lye, ammonia or even soda solution) in the alkaline range, may the pH-value be lowered respectively increased, so it counteracts the reaction of the sample (for example with chalk, kaolin, hydrated lime and clay). Suitable are here also Na₂HPO₄ (alkaline) or KH₂PO₄ (acidic) as a 0.1-1% solution.

2.2 Tips and Tricks

The larger the fine share of a sample, the greater is the dispersion effort. A necessary ultrasonic treatment lasting several minutes (or even longer) should be conducted in an external ultrasonic bath.

Here it is recommended to add to a 50 ml Erlenmeyer flask a spatula tip full (approximately 0.5-1 g) of sample and via simple wetting adding approximately 20 ml of the measuring liquid + the dissolved/ mixed dispersion aid (surfactant).

We recommend when utilizing surfactants those with low foam tendencies like for example Dusazin 901.

After briefly shaking the Erlenmeyer flask, it is fastened with laboratory clip inside the ultrasonic bath so the inner level of liquid is below the surface of the liquid of the ultrasonic bath.

After during the pre-tests determined deagglomeration period, the necessary suspension amount for the measurement will be added with a shaking motion via a pipette into the dispersion unit of the Laser Particle Sizer ANALYSETTE 22.

If the needed amount of solids exactly known, it may be weighed directly into the Erlenmeyer flask and after the dispersion the *complete* contents can be conveyed with a wash bottle, so no separation due to sampling with a pipette occurs.

With a few small tricks difficult samples like for example fly ash, sulphur, coal, plastics or pigments can be quickly dispersed even in water.

The sample is mixed with **one** drop of surfactant and after adding **one to two** drops of water with a glass stirrer or spatula mixed to a *paste*. Due to the relatively high share of surfactant compared to the share of water the cavitation is reduced and the sample quickly wetted.

After further addition of several drops of water and stirring at the same time it can be checked if not already wetted particles are swimming on top of the surface of the liquid.

Now the suspension is diluted down to approximately 20-30 ml and deagglomerated in the ultrasonic bath.

When using simple surfactants like dish soap it often occurs that when *pasting* and mixing the sample “foam” evolves which after the dispersion swims on top of the surface and is traversed into the particle sizer.

A too high stirring intensity pulls the foam into the measuring circuit and “*coarse particles*” are measured which are not present.

In order to avoid this “faulty measurement” it is possible to destroy the “*foam*” inside the Erlenmeyer flask: the tip of a glass stirrer is dipped into *n-butanol* so the glass stirrer is just wetted but no drops are recognizable.

By slightly touching the surface of the foam it collapses and the suspension can easily be analysed for particle size.

An additional possibility – with poor wetting – would be the addition of 2 to 3 drops of alcohol (for example ethyl alcohol) directly on the dry laboratory sample, which immediately absorbs the alcohol like a sponge. Now again water and a wetting agent can be added and dispersed accordingly.

The user, working with the “Small Volume Wet Dispersion Unit” of course has the possibility to utilize alcohol, alkanes, high-boiling benzines or other organic liquids.

It should be mentioned, that the Small Volume Wet Dispersion Unit is **not** explosion protected – please select the corresponding liquids – and only operate in well ventilated areas.

The exchange of the measurement liquid for example alcohols to other organic liquids is relatively easy, since many liquids can be mixed easily amongst each other.

The compatibility of solvents with the connective hoses should also be kept in mind. The seals in the measuring cell and the connective hoses are made of Viton.

Acetone (ketones), acetates and enamel thinner cannot be used.

The *resistance lists* available from hose manufacturers contain for the most common elastomers a rating of the chemical resistance against various operating mediums (liquids).

After a conducted measurement with the ANALYSETTE 22 it should become habit to *directly* rinse

the measuring system in order to avoid an unnecessary “depositing” or “adhering” of particles in the measuring circuit and especially on the measuring cell.

An „in-between rinsing“ with a surfactant is very helpful.

Residue on the measurement cell may be caused by tap water respectively limy water. Either it should be switched to distilled water or the measurement cell has to be cleaned from time to time.

The lime residue is removed in a few minutes by rinsing with 10% hydrochlorid acid. Afterwards it should be rinsed twice with regular water.

On the following pages you find a list of materials with the suitable dispersion liquids and additives. Please note the list is only available in German. Should you have questions please contact the FRITSCH laboratory: gerber@fritsch.de

Materials and suitable measuring liquids

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Author: Dipl. Chem. Ulrich Gerber, Fritsch GmbH, E-Mail: info@fritsch.de